

16mm width

19mm width

22mm width

25mm width

Size 16

Smaller diameter walking sticks and canes

Size 19

Most walking sticks

Size 22

Most adjustable crutches

Size 25

Most walking frames and shower chairs

Cut out the half-circles above, insert your cane, stick, crutches or frame to confirm the size you require.

australia flexyfoot size selector card



Note: Please measure shaft near the foot.

Please check that your printout is to scale by checking against a ruler

1cm

2cm

3cm

4cm

5cm

6cm

7cm

8cm

9cm

10cm

11cm

12cm

13cm

14cm

15cm